

MetaCell[®] Glucose solution (20%)

Product Description

MetaCell[®] Glucose Solution is a D-glucose solution with a concentration of 20% (w/w). Glucose is a key component that helps the growth and survival of mammalian cells under in vitro conditions and is crucial for cell yield and quality. MetaCell[®] Glucose Solution is used as a supplement for serum-free media suitable for cell culture.

This product is intended for research or further manufacturing but not for human or therapeutic use.

Product Name	Cat No.	Size	Storage	Shelf Life	Application
MetaCell [®] Glucose Solution (20%)	L5009-0100	100mL	2-8°C, protected from light	12 months	Mammalian cell culture

Product Applications

MetaCell[®] Glucose Solution can be added to a variety of MetaCell[®] cell culture media, including but not limited to MetaCell[®] DMEM, MetaCell[®] RPMI 1640, etc. CellPlus Bio offers multiple types of culture media to which MetaCell[®] Glucose Solution can be added.

Precautions

1. When using MetaCell[®] Glucose Solution, pay attention to aseptic operation to avoid contamination.
2. Each batch of MetaCell[®] Glucose Solution has undergone strict testing to ensure the consistency of product quality and performance between batches.